

# **Readme for dconv\_multiV4**

**A. Puri (puri@esrf.fr) 2 November 2023**

# Table of Contents

<b>1</b>	<b>dconv_multiV4 operation manual .....</b>	<b>1</b>
1.1	Intro .....	1
1.2	How to run the program .....	1
1.3	Column meaning .....	1
1.4	How to choose the right columns .....	1

# 1 dconv\_multiV4 operation manual

## 1.1 Intro

dconv\_multiV4 is an application for LISA data reduction. It writes an ASCII file as output containing the calibrated energy and the principal counters needed for data analysis.

## 1.2 How to run the program

Place the dconv\_configV4.txt configuration file in the same directory of the .exe file.

Double click the dconv\_multiV4.exe file, and wait for the program to load, it may take a while.

Select the file(s) to convert and press Convert.

The output files are written in the same directory of the input files and have the same name of the datafile with a final \_c.

## 1.3 Column meaning

eBraggEnergy: Calibrated energy in eV

I0\_EH2: Incoming x-ray intensity on the sample in EH2

I1\_EH2: Transmitted x-ray intensity after the sample in EH2

IX\_EH2: Transmitted x-ray intensity after the reference in EH2

IR\_EH2: Total Electron Yield counter

c8 and c9: free counters

fluo01, fluo02, . . . , fluo13: fluorescence counters

I0\_EH1: Incoming x-ray intensity on the sample in EH1

I1\_EH1: Transmitted x-ray intensity after the sample in EH1

IX\_EH1: Photodiode counter

## 1.4 How to choose the right columns

For transmission measure in EH1

absorption coefficient of the sample,  $\mu = \ln\left(\frac{I0\_EH1}{I1\_EH1}\right)$

For transmission measure in EH2

absorption coefficient of the sample,  $\mu = \ln\left(\frac{I0\_EH2}{I1\_EH2}\right)$

For fluorescence measure

normalized fluorescence,  $fluo\_norm = \left(\frac{fluo01+fluo02+...}{I0\_EH2}\right)$

Absorption coefficient of the reference,  $\mu\_ref = \ln\left(\frac{I1\_EH2}{IX\_EH2}\right)$